



- 源自芬兰经典设计
- 使用轻巧便捷，可减少手部疲劳
- 10µl以上单道配有疏水性滤芯
- 量程范围广(0.1µl-10000µl)
- 管嘴连件部分可高温灭菌

Technology Creates The Future

Genex 移液器

BIO-DL Genex 系列移液器

订购信息

产品编号	通道数	规格	可变量	测量体积	不准确性	不精确性	适用吸头
72441011	1	0.1-2.5µl	0.05µl	2.5µl	±2.50%	2.00%	10µl
				1.25µl	±3.00%	3.00%	
				0.25µl	±12.00%	6.00%	
72441031	1	0.5-10µl	0.1µl	10µl	±1.00%	0.80%	10µl
				5µl	±1.50%	1.50%	
				1µl	±2.50%	1.50%	
72441041	1	2-20µl	0.5µl	20µl	±0.90%	0.40%	200µl & 300µl
				10µl	±1.20%	1.00%	
				2µl	±3.00%	2.00%	
72441051	1	5-50µl	0.5µl	50µl	±0.60%	0.30%	200µl & 300µl
				25µl	±0.90%	0.60%	
				5µl	±2.00%	2.00%	
72441061	1	10-100µl	1µl	100µl	±0.80%	0.20%	200µl & 300µl
				50µl	±1.00%	0.40%	
				10µl	±3.00%	1.00%	
72441071	1	20-200µl	1µl	200µl	±0.60%	0.20%	200µl & 300µl
				100µl	±0.80%	0.30%	
				20µl	±2.50%	0.80%	
72441081	1	100-1000µl	5 µl	1000µl	±0.60%	0.20%	1000µl
				500µl	±0.70%	0.25%	
				100µl	±2.00%	0.70%	
72441101	1	1 - 5 ml	50 µl	5000µl	±0.50%	0.20%	5 ml
				2500µl	±0.60%	0.30%	
				1000µl	±0.70%	0.30%	
72441111	1	1-10ml	100 µl	10000µl	±0.60%	0.20%	10ml
				5000µl	±1.20%	0.30%	
				1000µl	±3.00%	0.60%	

产品编号	通道数	规格	可变量	测量体积	不准确性	不精确性	适用吸头
72471141	6	50-300µl	5µl	300µl	±0.70%	0.25%	300µl & 350µl
				150µl	±1.00%	0.5%	
				50µl	±1.50%	0.8%	
72411031	8	0.5-10µl	0.1µl	10µl	±1.50%	1.50%	10µl
				5µl	±2.50%	2.50%	
				1µl	±4.00%	4.00%	
72411051	8	5-50µl	0.5µl	50µl	±1.00%	0.50%	200µl & 300µl
				25µl	±1.50%	1.00%	
				5µl	±3.00%	2.00%	
72411061	8	10-100µl	0.5µl	100µl	±0.80%	0.20%	200µl & 300µl
				50µl	±1.00%	0.40%	
				10µl	±3.00%	1.00%	
72411141	8	50-300µl	5µl	300µl	±0.70%	0.25%	300µl & 350µl
				150µl	±1.00%	0.50%	
				50µl	±1.50%	0.80%	
72411081	8	100-1000µl	10µl	1000µl	±0.60%	0.20%	1000µl
				500µl	±0.70%	0.25%	
				100µl	±2.00%	0.70%	
72421031	12	0.5-10µl	0.1µl	10µl	±1.50%	1.50%	10µl
				5µl	±2.50%	2.50%	
				1µl	±4.00%	4.00%	
72421051	12	5-50µl	0.5µl	50µl	±1.00%	0.50%	200µl & 300µl
				25µl	±1.50%	1.00%	
				5µl	±3.00%	2.00%	
72421061	12	10-100µl	0.5µl	100µl	±0.80%	0.20%	200µl & 300µl
				50µl	±1.00%	0.40%	
				10µl	±3.00%	1.00%	
72421141	12	50-300µl	5µl	300µl	±0.70%	0.25%	300µl & 350µl
				150µl	±1.00%	0.50%	
				50µl	±1.50%	0.80%	